



Incus
SURGICAL LTD

THE VITAL LINK



Cranial Stabilization System


Welcome to Incus Surgical Limited

For nearly 30 years, Incus Surgical Limited have built and sustained an enviable reputation for the manufacture and supply of surgical instrumentation, covering every type of surgical speciality.

We have distinguished ourselves as a major leader in redefining the provision of high quality surgical instruments, and have become well known for the passion of our people and partners in providing unsurpassed levels to our customers internationally.

Our mission is to be the healthcare industry's first choice for all their surgical instrument requirements, and to deliver products of outstanding quality with great service, at an affordable price.


We are delighted to present our new Cranial Stabilization System — for prices and stock availability please email marie@incus-surgical.com





THE INCUS CRANIAL STABILIZATION SYSTEM

- 2 designs available - with Adult Horseshoe Headrest or with Skull Clamp
- Manufactured from Lightweight Aluminium with Stainless Steel fittings
- Skull Pins manufactured from either Stainless Steel or Titanium
- Movable in all directions and height adjustable
- Compatible with most existing systems
- Can be combined with the Incus Neuro Retractor System



THE INCUS CRANIAL STABILIZATION SYSTEM WITH SKULL CLAMP



NI37-200 Cranial Stabilization System with Skull Clamp, Aluminium
Consisting of 1 x **NI37-215**, 1 x **NI37-217** and 1 x **NI37-219**
(Skull Pins not included - to be ordered separately - 3 Required)



THE INCUS

CRANIAL STABILIZATION SYSTEM

WITH SKULL CLAMP — ACCESSORIES



NI37-219 Skull Clamp Only, Aluminium
(Skull Pins not included - to be ordered separately - 3 Required)



NI37-221 Skull Pin, Adult, Stainless, Steel Reusable

NI37-223 Skull Pin, Child, Stainless, Steel Reusable

NI37-225 Skull Pin, Adult, Titanium, Reusable

NI37-227 Skull Pin, Child, Titanium, Reusable



NI37-215 Base Unit and
Connecting Arm Only,
Aluminium and Stainless Steel



www.incus-surgical.com

THE INCUS CRANIAL STABILIZATION SYSTEM— ACCESSORIES



NI37-217 Swivel Adaptor Only,Aluminium



THE INCUS CRANIAL STABILIZATION SYSTEM WITH HEADREST



NI37-210 Cranial Stabilization System with Headrest, Aluminium
Consisting of 1 x **NI37-215**, 1 x **NI37-217**, 1 x **NI37-242** and 1 x **NI37-244**
(Skull Pins not included - to be ordered separately - 3 Required)



THE INCUS CRANIAL STABILIZATION SYSTEM WITH HEADREST — ACCESSORIES



NI37-242 Horseshoe Headrest only, Aluminium



NI37-244 Headrest Adult Gel-Pad only



THE INCUS

CRANIAL STABILIZATION SYSTEM




N137-260 Retractor System for Operational Hand Rest and Accessories, Aluminium
Complete, comprising: 2 x **N137-264**, 2 x **N137-266**, 2 x **N137-268** and 1 x **N137-270**

N137-264 Slightly Curved Arms Only, Pair - Left and Right, Aluminium

N137-266 Support for Curved Arms Only, Aluminium

N137-268 Connecting Clamp Only, Aluminium

N137-270 Connecting Arms, Pair - Left and Right, Aluminium



THE INCUS

CRANIAL STABILIZATION SYSTEM—

TABLE ACCESSORIES



NI37-252 Cranial Stabilization System, Crossbar, Aluminium
Width 64cm-70cm
(Requires 2 x **NI37-254** Universal Slide Rail Fittings)



NI37-254 Cranial Stabilization System, Universal Slide Rail Fitting




www.incus-surgical.com

THE INCUS CRANIAL STABILIZATION SYSTEM— TABLE ACCESSORIES



NI37-250 Table Rail Adaptor, Aluminium



Contact:

Incus Surgical Limited

Mabrook House

Bocking End

Braintree

Essex

CM7 9AA

United Kingdom

T: +44 (0) 1376 347300

E: marie@incus-surgical.com

www.incus-surgical.com